

Press Release

NEOFORCE LAUNCHES NeoPEEP™ NEONATAL RESUSCITATION CIRCUIT

IVYLAND, PA USA:, THURSDAY, AUGUST 16, 2007.

NeoForce is proud to introduce the NeoPEEP™ Neonatal Resuscitation Circuit. “Finally we have a disposable single patient use circuit to meet all of your resuscitation needs” stated Otho Boone, President of NeoForce Group. The NeoPEEP™ Circuit features Positive End Expiratory Pressure control (PEEP) which allows you to manage Functional Residual Capacity (FRC). The circuit will provide consistent PEEP throughout the resuscitation process and is designed for use on patients weighing less than 10kg.

Our NeoPeep™ circuit includes a T-piece connector with a peep control cap allowing adjustments from 0-45cm of water. Inspiration and expiration can be easily controlled with the finger/thumb hole on the end of the cap. The circuit hose has a length of 59 inches with a standard 15mm O.D. device connector. The additional length will come in handy stated Boone, “many nurseries are very space constrained and the extra length will give users more flexibility to reach both the resuscitation equipment and the patient.” NeoForce also includes a “Soft Adapter” with each circuit to allow quick connection to other resuscitation devices. The NeoPeep™ Circuit is available with pink or blue T-Pieces and is sold in quantities of twenty circuits per box.

Contact: Nancy St Clair VP Marketing
Phone: (215) 672-6800 x 202
e-mail: nancy.stclair@neoforcegroup.com

FOR IMMEDIATE RELEASE

-End-